

# CT Tibia-Fibula without contrast

Clinical Indication	Comment
Trauma, pre-operative planning, pain	

## ACQUISITION

<b>Contrast Phase</b>	
<b>Extent</b>	
<b>Scan direction</b>	
<b>Acquisition technique</b>	
<b>IV Contrast technique</b>	NA
<b>IV Contrast rate</b>	NA
<b>Saline flush</b>	
<b>Delay</b>	
<b>SmartPrep/ Bolus Pro</b>	
<b>Patient Age</b>	All ages
<b>Patient Size</b>	All sizes
<b>Dose Index Range</b>	Action level for investigation: >500 mGy total CTDIvol, >5000 mGy*cm total DLP

## RECONSTRUCTION

Source	Name of series	Reconstruction/Reformat
Raw data	SAG	Sagittal
Raw data	COR	Coronal
Raw data	3D	3D Volume rendered

## Keywords

Trauma  
pre-operative planning  
pain